# CALPROTECTIN FOR LEDDVÆSKE

**HURTIGTEST FOR RISIKOSTRATIFIKASJON AV  
INFEKSJON I PASIENTER MED LEDDPROTESE**

**BRUKSANVISNING**

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Personal information   Clinic information   Summary

**Contact person for the clinic**  
This person will serve as a Clinic Administrator

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User types in Lyfstone® CLP system	Clinic Admin	Lab Tester	Bedside Tester
Register Clinics	X		
Register Lab Tester	X		
Register Bedside Tester	X		
View all patient's test results	X	X	X
Web access	X	X	X
Perform test on Lyfstone® app	X	X	X

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### Create staff user

First name* <small>Required</small>	Language* <small>Required</small>
Last Name <small>Required</small>	Clinic* <small>Required</small>
Email address* <small>Required</small>	<input checked="" type="checkbox"/> Is Active
Country	Roles
Address line 1	<input type="checkbox"/> CLINIC_ADMINISTRATOR
Address line 2	<input type="checkbox"/> CLINIC_LAB
Primary phone    City	<input type="checkbox"/> BEDSIDE
Secondary phone	

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The screenshot shows a table titled "Manage results" with a search bar and a "SEARCH" button. The table has columns for Date, System id, Sample id, Clinic, Batch, Concentration, Classification, and DETAILS. The Classification column contains red "CLINIC" labels and an orange "BEDSIDE" label. The bottom right of the table shows "Rows per page: 25" and "1/3 of 3".

Date	System id	Sample id	Clinic	Batch	Concentration	Classification	DETAILS
11.5.2022, 11:02:01	000001	None	Kennel Clinic 2	20222009	564	CLINIC	DETAILS
11.5.2022, 11:04:56	000002	None	Kennel Clinic 2	20222009	531	CLINIC	DETAILS
13.5.2022, 09:41:34	000005	001	Kennel Clinic 2	20222009	474	BEDSIDE	DETAILS

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## 16. REFERANSER:

1. **G Hetland, G J Talgö, M K Fagerhol** *Chemotaxins C5a and fMLP induce release of calprotectin (leucocyte L1 protein) from polymorphonuclear cells in vitro* J Clin Pathol: Mol Pathol 1998;51:143–148
2. **Cuida, M. Brun, J.G., Johannessen, A. C., Jonsson, R. (1996)** *Immunohistochemical characterization of the cellular infiltrates in Sjögren's syndrome, rheumatoid arthritis and osteoarthritis with special reference to calprotectin-producing cells* APMIS, 104, 7-8, 881-890  
<http://dx.doi.org/10.1111/j.1699-0463.1996.tb04953.x>
3. **Cesaro A, Anceriz N, Plante A, Pagé N, Tardif MR, Tessier PA (2012)** *An Inflammation Loop Orchestrated by S100A9 and Calprotectin Is Critical for Development of Arthritis.* PLoS ONE 7(9): e45478. <https://doi.org/10.1371/journal.pone.0045478>
4. **Teng, T.-S., Ji, A., Ji, X.-Y., & Li, Y.-Z. (2017).** *Neutrophils and Immunity: From Bactericidal Action to Being Conquered.* Journal of Immunology Research, 2017, 9671604. <http://doi.org/10.1155/2017/9671604>
5. **STRÍŽ, I., TREBICHAŤSKÝ, I., (2004)** *Calprotectin- a pleiotropic Molecule in Acute and Chronic inflammation.* Physiol. Res. 53:245-253
6. **Urban CF, Ermert D, Schmid M, Abu-Abed U, Goosmann C, Nacken W, et al. (2009)** *Neutrophil Extracellular Traps Contain Calprotectin, a Cytosolic Protein Complex Involved in Host Defense against Candida albicans.* PLoS Pathog 5(10): e1000639. <https://doi.org/10.1371/journal.ppat.1000639>
7. **M. Wouthuyzen-Bakker, J. J. W. Ploegmakers, G. A. Kampinga, L. Wagenmakers-Huizenga, P. C. Jutte, A. C. Muller Kobold** *Synovial calprotectin: a potential biomarker to exclude a prosthetic joint infection*  
DOI: <http://10.1302/0301-620X.99B5.BJJ-2016-0913.R2> Published 28 April 2017
8. **C. Nich, Y. Takakubo, J. Pajarinen, M. Ainola, A. Salem, T. Sillat, A.J. Rao, M. Raska, Y. Tamaki, M. Takagi, Y.T. Kontinen, S.B. Goodman, J. Gallo** *Macrophages—key cells in the response to wear debris from joint replacements* J. Biomed. Mater. Res. A, 101 (2013), pp. 3033-3045
9. **Bitar, D., & Parvizi, J. (2015).** *Biological response to prosthetic debris.* World Journal of Orthopedics, 6(2), 172–189. <http://doi.org/10.5312/wjo.v6.i2.172>
10. **Jared A Warren 1, Alison K Klika 1, Kathleen Bowers 1, Jessica Colon-Franco 1, Nicolas S Piuze 1, Carlos A Higuera 2** *Calprotectin Lateral Flow Test: Consistent Across Criteria for Ruling Out Periprosthetic Joint Infection.* J Arthroplasty. 2022 Feb 3; S0883-5403(22)00099-7. doi:10.1016/j.arth.2022.01.082.
11. **Lazic I, Prodinger P, Stephan M, Haug AT, Pohlif G, Langer S, von Eisenhart-Rothe R, Suren C.** *Synovial calprotectin is a reliable biomarker for periprosthetic joint infections in acute-phase inflammation - a prospective cohort study.* Int Orthop. 2022 May 7. doi: 10.1007/s00264-022-05421-1. Epub ahead of print. PMID: 35524793.

## 17. BESTILLINGSINFORMASJON

Product code: LYFCLP005

Lyfstone® Calprotectin for Synovial Fluid

Order at you local distributor or at [sales@lyfstone.com](mailto:sales@lyfstone.com)

Produsert av

CALPRO AS



Arnstein Arnebergsvei 30

N-1366 Lysaker, Norway

Net: [www.calpro.no](http://www.calpro.no)

Email: [mail@calpro.no](mailto:mail@calpro.no)

Phone: +47 40 00 42 79

<b>Symbols Key / Symbolschlüssel / Explication des symboles</b>	
	In Vitro Diagnostic Medical Device / In Vitro Diagnosticum / Dispositif médical de diagnostic in vitro
	Lot Number / Chargenbezeichnung / Numéro de lot
	Expiration Date / Verfallsdatum / Date de péremption
	Storage Temperature / Lagertemperatur / Température de conservation
	CE Mark / CE-Zeichen / Marquage CE
	Catalogue Number / Katalog Nummer / Référence du catalogue
	Contains sufficient for "n" tests / Ausreichend für "n" Tests / Contenu suffisant pour "n" tests
	Do not reuse / Nicht wiederverwenden / No reutilizar
	Manufactured by / Hergestellt von / Fabricado por
	Follow operating instructions / Bedienungsanweisung beachten / Siga las instrucciones de funcionamiento

Produsert av:

**CALPRO AS**

**Arnstein Arnebergs vei 30,**

**N-1366 Lysaker, Norway**

Web: [www.calpro.no](http://www.calpro.no)

Epost: [mail@calpro.no](mailto:mail@calpro.no)

Telefon: **+47 40 00 42 79**